


Banded Bypass with Ring




Banded Sleeve with Ring



03-ENG-2012


# Stable Weight loss in Banded Bypass and Banded Sleeve MiniMizer Ring



%EWL for standard RYGBP and its two hybrid forms of LL-RYGBP and Banded RYGBP (\* P<0.05).

Obrien PE, McPhail T, Chaston TB et al. Systematic review of medium-term weight loss after bariatric operations. *Obes Surg* 2006;16:1032-40

## MiniMIZER RING



Closing positions from the largest to the smallest ring size:

- 1 size 1 approx. 8.0 cm length, approx. 26mm internal diameter
- 2 size 2 approx. 7.5 cm length, approx. 24mm internal diameter
- 3 size 3 approx. 7.0 cm length, approx. 22mm internal diameter
- 4 size 4 approx. 6.5 cm length, approx. 20mm internal diameter



46 Albert Street  
 Preston VIC 3072  
 Phone 1300 244 473  
 Fax 03 8470 5490

AUS: 1300 244 473  
 NZ: 0800 424 769  
 www.biggreen.com.au  
 customerservice@biggreen.com.au



The difference makes it happen.  
MiniMizer Ring




- Placement
- 4-6 cm below the gastro-esophageal-junction
  - or 2 cm above the anastomosis



- The MiniMizer Ring is radiopaque
- needle and ring are well visible under X-ray
  - simple placement and closing
  - two loops for fixation after placement
  - can be reopened and closed in different diameter



MiniMIZER RING



Why Banded Bypass and Banded Sleeve?

- protection against pouch dilatation
- protects the anastomosis and the jejunum against overstretching
- minimizes dumping syndrome
- stabilizes long term weight loss

