



# Vathin<sup>®</sup> Single-use Flexible Video Rhinolaryngoscope

# Vathin<sup>®</sup> Single-use Flexible Video Rhinolaryngoscope





## Slim Diameter

A slim insertion tube allows easier approach to the nasal, pharyngeal and laryngeal cavities, making the scope of observation more comprehensive.

## Field of View

The field can be thoroughly inspected with the scope's wide 110° field of view.



## Bending Range

A  $\pm 210^\circ$  wide angulation range allows a smoothly insertion into the upper airway and offers the flexibility to perform better during the procedure.

## Specifications

Model	Color	Working channel ID(mm)	Insertion tube OD(mm)	Working length(mm)	Bending range	Field of view
RL-E1800	Grey	0	3.2	300	Up 0°~210° Down 0°~210°	110°± 15%
RL-E1E00	Green	2.2	4.9			



1st floor, Building 12 , No.9, West Chuanqi Road, Jiu Hua Economic and Technological Development Zone , Xiangtan 411100, China  
 Tel: +86-021-34781226 | Mail: [marketing@vathin.com](mailto:marketing@vathin.com) | [www.vathin.com](http://www.vathin.com) | Version\_MD22013



Distributed by Big Green:  
 1300 244 473 (AU) 0800 424 769 (NZ) | [customerservice@biggreen.com.au](mailto:customerservice@biggreen.com.au) | [biggreen.com.au](http://biggreen.com.au)