



Vathin[®] Single-use Flexible Video Rhinolaryngoscope

Vathin[®] Single-use Flexible Video Rhinolaryngoscope





Slim Diameter

A slim insertion tube allows easier approach to the nasal, pharyngeal and laryngeal cavities, making the scope of observation more comprehensive.



Field of View

The field can be thoroughly inspected with the scope's wide 110° field of view.



Bending Range

A $\pm 210^\circ$ wide angulation range allows a smoothly insertion into the upper airway and offers the flexibility to perform better during the procedure.

Specifications

Model	Color	Working channel ID(mm)	Insertion tube OD(mm)	Working length(mm)	Bending range	Field of view
RL-E1800	Grey	0	3.2	300	Up 0°~210° Down 0°~210°	110° \pm 15%
RL-E1E00	Green	2.2	4.9			



1st floor, Building 12, No.9, West Chuanqi Road, Jiu Hua Economic and Technological Development Zone, Xiangtan 411100, China
Tel: +86-021-34781226 | Mail: marketing@vathin.com | www.vathin.com | Version_MD22013



Distributed by Big Green:

1300 244 473 (AU) 0800 424 769 (NZ) | customerservice@biggreen.com.au | biggreen.com.au