



# PinkProtect™ Kit and Components

## Patient Positioning Solutions for the ICU and Critical Care Unit

Begin with the base kit and order additional components to address the positioning requirements of your patients. See the reverse side for details.

### Patient Safety Advantages

- Proven, skin-friendly positioners designed to address the specific patient safety challenges of the Intensive Care and Critical Care Units
- Featuring DermaProx Pink™ – a breathable, soft, moisture wicking material for the utmost protection against HAPI formation
- Supportive, high-density cores that help distribute pressure during lengthy patient positioning intervals
- Unique designs tailored to the specific requirements of Intensive Care and Critical Care Units


#### SideProne Headrest

Integrated ET Tube channel to maintain unobstructed airways; materials optimized for both skin protection and long-term support



#### Dual Layer Pads

Multi-purpose pads protect high-risk pressure areas with DermaProx Pink™ over a high density layer that enhances pressure distribution and durability



#### Pink Arcs

Maintain skin integrity on heels or tops of feet throughout lengthy positioning intervals.



Refer to hospital protocol when positioning patient and selecting positioning products.



# PinkProtect™ Kit and Components

Patient Positioning Solutions for the ICU and Critical Care Unit




**XD40544**

Pink Arcs (standard)  
16/case  
22.9cm x 20.3cm x 10.2cm




**XD40524**

Pink Arcs (medium)  
8/case  
40.6cm x 15.2cm x 10.2cm




**XD40545**

Pink Arcs (large)  
6/case  
40.6cm x 15.2cm x 15.2cm




**XD40770**

SideProne Headrest 12/  
case  
25.4cm x 26.6cm x 6.3cm



**XD40521**

DermaProx Pink  
12/case  
24.1cm x 26.6cm x 16.5cm




**XD40771**

Large Dual Layer Pad 24/  
case  
38.1cm x 27.9cm x 2.54cm



**XD40772**

Standard Dual Layer Pad 48/case  
- 8" x 11" x 1"



**XD40780**

Pink Protect Kit includes:  
XD40544, XD40770, XD40771,  
XD40772. 4/case