

The Tube takes you to Pylorus Station

Calibration Tube for Sleeve Gastrectomy, Gastric Bypass and Mini Gastric Bypass



TUBE 32

TUBE 34

TUBE 36

TUBE 38

TUBE 40



☎ AU 1300 244 473

☎ NZ 0800 424 769

📍 46 Albert Street Preston VIC 3072 Australia


🌐 www.biggreen.com.au

✉ customerservice@biggreen.com.au

Temperature stable silicone tube

Technical Information:

- 1 soft conical radiopaque tip with suction hole
- 2 6 small holes avoid tissue suction damage
- 3 cm marking for precise placement
- 4 universal suction connector
- 5 clip to close tube for leak test
- 6 colour coded



TUBE32W	32 Fr
TUBE34R	34 Fr
TUBE36B	36 Fr
TUBE38G	38 Fr
TUBE40Y	40 Fr

