

The Tube takes you to Pylorus Station

Calibration Tube for Sleeve Gastrectomy, Gastric Bypass and Mini Gastric Bypass



TUBE 32

TUBE 34

TUBE 36

TUBE 38

TUBE 40



☎ AU 1300 244 473

☎ NZ 0800 424 769

📍 46 Albert Street Preston VIC 3072 Australia


🌐 www.biggreen.com.au

✉ customerservice@biggreen.com.au

Temperature stable silicone tube

Technical Information:

- ① soft conical radiopaque tip with suction hole
- ② 6 small holes avoid tissue suction damage
- ③ cm marking for precise placement
- ④ universal suction connector
- ⑤ clip to close tube for leak test
- ⑥ colour coded



TUBE32W 32 Fr
TUBE34R 34 Fr
TUBE36B 36 Fr
TUBE38G 38 Fr
TUBE40Y 40 Fr

